CANADA USE ONLY

100481072.3

Instructions for Use

Cordis VISTA BRITE TIP[®] Guiding Catheters

Mode d'emploi

Cathéters-guide VISTA BRITE TIP[®] Cordis

Explanation of symbols on labels and packaging - CANADA USE ONLY

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Cordis VISTA BRITE TIP® Guiding Catheters

CANADA USE ONLY

STERILE. Sterilized with ethylene oxide gas. Nonpyrogenic. Radiopaque. For one use only. Do not resterilize.

Device Description

The Cordis Guiding Catheter has a large nontapered lumen that allows contrast medium injections and facilitates the intravascular passage of interventional devices.

Indications

The guiding catheter is intended for use for intravascular introduction of interventional/diagnostic devices into the coronary or peripheral vascular systems.

Contraindications

None known for guiding catheters.

Warnings

This product is designed and intended for single use. It is not designed to undergo reprocessing and re-sterilization after initial use. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.

Do not use with Ethiodol or Lipiodol* contrast media, or other such contrast media which incorporates the components of these agents.

Precautions

- Store in a cool, dark, dry place.
- Do not use open or damaged packages.
- Use prior to the "Use By" date.
- Do not resterilize.
- Do not expose to organic solvents.
- Inspect the guiding catheter before use to verify that its size, shape and condition are suitable for the specific procedure.
- If strong resistance is met during manipulation, discontinue the procedure and determine the cause of the resistance before proceeding. If the cause of the resistance can not be determined, withdraw the catheter.
- Torquing the guiding catheter excessively while kinked may cause damage which could result in possible separation along the catheter shaft. Should the guiding catheter shaft become severely kinked, withdraw the entire system (guiding catheter, guidewire and catheter sheath introducer).
- Advancement, manipulation and withdrawal of the guiding catheter should always be performed under fluoroscopic guidance.
- Extreme care must be taken to avoid damage to the vasculature through which the guiding catheter passes. The guiding catheter may occlude smaller vessels. Care must be taken to avoid complete blood flow blockage.
- Large internal lumen guiding catheters require less force on the syringe during injection.

Complications

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications.

Complications may occur at anytime during or after the procedure.

Possible complications include, but are not limited to the following:

- air embolism
- hematoma at the puncture site
- infection
- perforation of the heart or vessel wall

Recommended Procedure

1. Remove the guiding catheter from its packaging.
2. Inspect the guiding catheter upon removal from packaging to verify that it is undamaged.
Warning: Do not use a guiding catheter that has been damaged in any way. If damage is detected, replace with an undamaged guiding catheter.
3. Prior to use, flush the guiding catheter lumen with a heparinized saline solution.
4. Introduce the guiding catheter into the vasculature through the catheter sheath introducer, and/or over an indwelling guidewire using a percutaneous entry technique of choice.
5. Set up a continuous heparinized saline flush through the sidearm of a hemostatic valve attached to the guiding catheter hub.
Note: It is recommended that a continuous heparinized saline flush be maintained between the guiding catheter and any intraluminal device passed through it.
6. Under fluoroscopic guidance, advance the guiding catheter over the guidewire or introducer until desired position is attained.
7. Remove the guidewire or introducer prior to the introduction of other intravascular devices or infusion of contrast agents.

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

THERE IS NO EXPRESS OR IMPLIED WARRANTY, INCLUDING WITHOUT LIMITATION ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ON THE CORDIS PRODUCT(S) DESCRIBED IN THIS PUBLICATION. UNDER NO CIRCUMSTANCES SHALL CORDIS BE LIABLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES OTHER THAN AS EXPRESSLY PROVIDED BY SPECIFIC LAW. NO PERSON HAS THE AUTHORITY TO BIND CORDIS TO ANY REPRESENTATION OR WARRANTY EXCEPT AS SPECIFICALLY SET FORTH HEREIN.

Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Cordis Corporation will not be responsible for any direct, incidental or consequential damages resulting from reuse of the catheter.

Protected under the following U.S. Patents: 5,045,072 and 5,171,232.

* Ethiodol and Lipiodol are trademarks of Guerbet S.A.

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© 2018, Cardinal Health. All Rights Reserved. CORDIS, the Cordis LOGO, CARDINAL HEALTH and VISTA BRITE TIP are trademarks of Cardinal Health and may be registered in the US and/or in other countries.
* Ethiodol et Lipiodol sont des marques de commerce de Guerbet S.A.